

# PRODUCT SPECIFICATION

## SAFETWARE HDPE Glove

- Cost-effective choice among disposable gloves, providing an economical solution for various tasks.
- Odorless and food-safe, making them suitable for handling food items without affecting taste or quality.
- Ideal for tasks requiring minimal barrier protection, offering basic protection for low-risk activities.
- Lightweight design ensures quick and easy donning, enhancing user comfort and efficiency.
- Allergy-free formulation, completely free from Type I Latex and Type IV allergies, reducing the risk of allergic reactions.
- Disposable and hygienic, designed for single use to prevent cross-contamination.
- Ambidextrous fit, easily worn on either hand for convenience.

## STORAGE

- Store in a cool, dry place away from direct sunlight or heat sources to prevent material degradation.
- Keep in original packaging until ready for use to maintain cleanliness and hygiene.
- Store at room temperature, ideally between 10°C to 40°C, to preserve glove flexibility and durability.
- Check for damage before use, such as tears or holes, to ensure effective protection.

## TECHNICAL DATA

MATERIAL	High Density Polyethylene (HDPE)
COLOUR	Clear
SIZE	Universal
LENGTH	270 ± 10 mm
WEIGHT (PER PIECE)	0.7 gram
PACKAGING	100 pcs/ Printed polybag, 100 bags/ carton
STANDARD	EN 455 part 2, EN 455-1, EN 374-2
SHELF LIFE	5 Years

Item Code: PEG07



## APPLICATIONS

- Cleaning tasks
- Food preparation and handling
- Pet care
- Arts and crafts
- Hair dyeing and salon use