

# PRODUCT SPECIFICATION

## SAFETYWARE NEX Powder Free Nitrile Exam Gloves 3.5g Finger Textured

- Medical grade of AQL 1.5, suitable for healthcare and dental applications.
- It is latex-free, which could eliminate latex-allergy concerns.
- It is lightweight for enhanced dexterity and comfort.
- Textured grip for secure handling and precision.
- Conforms to ASTM D6319, ASTM D412.

### APPLICATIONS

- Laboratory
- Cleaning Work
- General Industry

### STORAGE

- Store in a cool, dry place. Avoid excessive heat.
- Open box should be shielded from exposure to direct sun or fluorescent lighting.

### TECHNICAL DATA

MATERIAL	100% Synthetic Nitrile	
SURFACE TEXTURE	Beaded Cuff, Finger Textured	
COLOR	Blue	
SIZE	S, M, L, XL	
WEIGHT	3.5 ± 0.3 g	
THICKNESS (SINGLE WALL)	Finger : 0.10 ± 0.02 mm Palm : 0.08 ± 0.02 mm	
	Before Aging	After Aging
TENSILE STRENGTH (Mpa)	Min 14	Min 14
ELONGATION AT BREAK (%)	500	400
LENGTH	240 ± 5 mm	
STANDARD	ASTM D6319, ASTM D412	
PACKAGING	100 pcs /box, 10 boxes /carton	
SHELF LIFE	3 Years	

Item Code: NEX35



### ITEM CODE

SIZE	ITEM CODE	PALM WIDTH (± 5.0 mm)
S	NEX35-BLU-SM	85
M	NEX35-BLU-MD	95
L	NEX35-BLU-LG	105
XL	NEX35-BLU-XL	115

**SAFETYWARE SDN. BHD.** (632971-X)

Plot 237, Lengkok Perindustrian Bukit Minyak 3,  
Bukit Minyak Industrial Estate,

14100 Simpang Ampat, Penang, Malaysia

Tel : +604-5023 882 | Fax : +604-5081 882  
www.safetyware.com | info@safetyware.com

It is the product users' responsibility to conduct all the appropriate evaluations & testing to determine the suitability of SAFETYWARE products for use in a particular environment or purpose. SAFETYWARE disclaims all warranties other than as expressly provided.

Date created: 14/11/2023

Revision: 02

Date revised: 10/07/2024